


Product Datasheet

Mouse M-CSF R / CSF1R / CD115 Protein (orb383507)

Description

Mouse M-CSF R, His Tag is expressed from human



SDS-PAGE
analysis of
Mouse M-
CSF protein

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

M-CSF R

Preservatives


PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.



Immobilized
Mouse M-
CSF Protein
at 5 µg...

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.1 kDa. The protein migrates as 66-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Mouse M-CSF R, His Tag (orb383507) is expressed from human 293 cells (HEK293). It contains AA Ala 20 - Ser 511 (Accession # NP_001032948).

MW

57.1 kDa

Uniprot ID

[P09581](#)

NCBI

[NP_001032948](#)

Expiration Date

6 months from date of receipt.