

Product Datasheet

Mouse PD-L2 / B7-DC Protein (orb383511)

Description

Mouse PD-L2, Fc Tag (PD2-M5254) is expressed from


Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

PD-L2

Preservatives

Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 49.2 kDa. As a result of glycosylation, the protein migrates as 60-70 kDa under reducing (R) condition, and 120-140 kDa under non-reducing (NR) condition (SDS-PAGE).

Protein Sequence

NP_067371

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Mouse PD-L2, Fc Tag (orb383511) is expressed from human 293 cells (HEK293). It contains AA Leu 20 - Arg 219 (Accession # NP_067371).

MW

49.2 kDa

Uniprot ID

[Q9WUL5](#)

NCBI

[NP_067371](#)

Expiration Date

6 months from date of receipt.

SDS-PAGE
analysis of
Mouse PD-
L2 protein