

Product Datasheet

ENaC alpha Antibody (orb383746)

Description

Mouse Anti-Rat ENaC alpha Monoclonal IgG2b

Species/Host

Mouse

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

IHC, WB

Applications

Immunogen

Synthetic peptide from the N-terminal of Rat ENaC alpha (aa. 46-68)

Target

ENaC alpha

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

-20°C

Note

For research use only

Application notes

A 1:1000 dilution of SMC-239 was sufficient for detection of ENaC alpha in 15 µg of Mouse whole kidney lysate by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2b

Clonality

Monoclonal

MW

85kDa

Uniprot ID

[Q6IRJ1](#)

NCBI

[NP_113736](#)

Entrez



[25122](#)

Dilution Range


WB (1:1000); IHC-P (1:150)

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of mouse
Ki...



Western Blot analysis of mouse whole kid...