

Product Datasheet

ENaC beta antibody (orb383768)

Description

Mouse monoclonal antibody against ENaC

Species/Host	Mouse
Reactivity	Mouse
Conjugation	ATTO-594
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the C-terminal of Rat ENaC beta (aa. 617-638)
Target	ENaC beta
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-240 was sufficient for detection of ENaC beta in 15 µg of Mouse whole kidney lysate by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	87kDa
Uniprot ID	P37090
NCBI	NP_036780
Entrez	24767
Dilution Range	WB (1:1000); IHC (1:150)
Expiration Date	12 months from date of receipt.