


Product Datasheet

ENaC alpha Antibody (orb383800)

Description	Mouse Anti-Rat ENaC alpha Monoclonal IgG1
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide from the N-terminal of Rat ENaC alpha (aa. 46-68)
Target	ENaC alpha
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-242 was sufficient for detection of ENaC alpha in 15 µg of Mouse whole kidney lysate by ECL immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	85kDa
Uniprot ID	Q6IRJ1
NCBI	NP_113736
Entrez	25122
Dilution Range	WB (1:1000); IHC-P (1:150)
Expiration Date	12 months from date of receipt.



Western Blot analysis of mouse whole kidney...