

Product Datasheet

ATG4B antibody (orb383937)

Description	Rabbit polyclonal antibody against ATG4B conjugated to Biotin.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the N-terminal to mid of Human ATG4B (aa. 94-106)
Target	ATG4B
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-630 was sufficient for detection of ATG4B in 15 µg of Human HeLa Cell Lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	44 kDa
Uniprot ID	Q9Y4P1
NCBI	NP_037457.3
Entrez	23192
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.