



Product Datasheet

RTP3 antibody (orb38411)


Description	RTP3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Receptor-transporting protein 3 protein (1-211AA)
Target	RTP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BQQ7
Dilution Range	WB:1:2000-8000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Mouse spleen t...



, Mouse lung tissue, Mouse kidney tissue...