

Product Datasheet

HCK antibody (orb38421)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	HCK
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59600
Uniprot ID	P08631
NCBI	NP_001165604.1 , NP_001165600.1 , NP_002101.2 , NP_001165602.1 , NP_001165601.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

