

Product Datasheet

Calcium Sensing Receptor antibody (orb384333)

Description

Rabbit polyclonal antibody against CASR

Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ICC, WB
Immunogen	Synthetic peptide from the N-terminal (1-100 aa) of Human calcium-sensing receptor (CASR)
Target	Calcium Sensing Receptor
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-713 was sufficient for detection of calcium-sensing receptor (CASR) in 15 µg of Human HeLa Cell Lysates by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	140 kDa
Uniprot ID	P41180
NCBI	NP_000379.2
Entrez	846
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.