

Product Datasheet

HMGB1 protein (orb384348)

Description	HMGB1 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	HMGB1
Form/Appearance	0.14mol/L NaCl, 0.027mol/L KCl, 0.01mol/L Na ₂ HPO ₄ ,0.018mol/L KH ₂ PO ₄ ,pH7.4
Concentration	1mg/ml
Note	For research use only
Application notes	N-Trx-his6
Purity	> 90%
MW	38.5 kDa
Uniprot ID	P09429
NCBI	002128
Expiration Date	6 months from date of receipt.