

## Product Datasheet

### HPV16E6 protein (orb384350)

|                            |                                |
|----------------------------|--------------------------------|
| <b>Description</b>         | HPV16E6 protein                |
| <b>Reactivity</b>          | Virus                          |
| <b>Conjugation</b>         | Unconjugated                   |
| <b>Tested Applications</b> | ELISA, SDS-PAGE, WB            |
| <b>Target</b>              | HPV16E6                        |
| <b>Form/Appearance</b>     | 10mM Tris, pH 8.0              |
| <b>Concentration</b>       | 0.5 mg/ml                      |
| <b>Note</b>                | For research use only          |
| <b>Application notes</b>   | pET28a                         |
| <b>Purity</b>              | > 85%                          |
| <b>MW</b>                  | 23.5 kDa                       |
| <b>Uniprot ID</b>          | <b>Q77E16</b>                  |
| <b>Expiration Date</b>     | 6 months from date of receipt. |