

Product Datasheet

CD304 protein (orb384351)

Description	CD304 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	CD304
Form/Appearance	50mM Tris, 50mM NaCl, pH 8.0
Concentration	0.5 mg/ml
Note	For research use only
Application notes	N-His6
Purity	> 85%
MW	21.6 kDa
Uniprot ID	O14786
NCBI	007533
Expiration Date	6 months from date of receipt.