


Product Datasheet

SPZ1 antibody (HRP) (orb38483)

Description	SPZ1 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Spermatogenic leucine zipper protein 1 protein (181-430AA)
Target	SPZ1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXG8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse gonadal t...