

Product Datasheet

NTF3 antibody (orb38486)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	NTF3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	29355
Uniprot ID	P20783
NCBI	NP_002518.1 , NP_001096124.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

