


Product Datasheet

Human ACAN ELISA Kit (orb385264)

Description	Sandwich ELISA for quantitative detection of human Aggrecan concentrations.
Reactivity	Human
Range	1560 pg/ml - 100000 pg/ml
Tested Applications	ELISA
Target	ACAN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Urine
Uniprot ID	P16112
Sensitivity	50 pg/ml
Entrez	176
Expiration Date	Please enquire.



Representative
Standard
Curve