

Product Datasheet

Human ADAM8 ELISA Kit (orb385268)

Description Sandwich ELISA for quantitative detection of human CD156a concentrations.

Reactivity Human

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target ADAM8

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum

Uniprot ID **P78325**

Sensitivity 10 pg/ml

Entrez **101**

Expiration Date Please enquire.



Representative
Standard
Curve