

Product Datasheet

Human CD97 ELISA Kit (orb385270)

Description Sandwich ELISA for quantitative detection of human CD97 concentrations.

Reactivity Human

Range 125 pg/ml - 8000 pg/ml

Tested Applications ELISA

Target CD97

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P48960**

Sensitivity 10 pg/ml

Entrez **976**

Expiration Date Please enquire.



Representative
Standard
Curve