


Product Datasheet

Human Apolipoprotein E ELISA Kit (orb385274)

Description	Human Apolipoprotein E ELISA Kit
Reactivity	Human
Range	3120 pg/ml - 200000 pg/ml
Tested Applications	ELISA
Target	APOE
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Saliva, Urine, Human Milk
Uniprot ID	P02649
Sensitivity	50 pg/ml
Entrez	348
Expiration Date	Please enquire.



Representative
Standard
Curve