


Product Datasheet

Human Caspase 8 ELISA Kit (orb385277)

Description	Sandwich ELISA for quantitative detection of human Caspase 8 concentrations.
Reactivity	Human
Range	78 pg/ml - 5000 pg/ml
Tested Applications	ELISA
Target	CASP8
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Tissue Homogenates
Uniprot ID	Q14790
Sensitivity	35 pg/ml
Entrez	841
Expiration Date	Please enquire.



Representative
Standard
Curve