

Product Datasheet

Human CD44 ELISA Kit (orb385281)

Description Sandwich ELISA for quantitative detection of human CD44 concentrations.

Reactivity Human

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target CD44

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P16070**

Sensitivity 80 pg/ml

Entrez **960**

Expiration Date Please enquire.



Representative
Standard
Curve