


Product Datasheet

Human EDIL3 ELISA Kit (orb385290)

Description	Sandwich ELISA for quantitative detection of human EDIL3 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	EDIL3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	O43854
Sensitivity	10 pg/ml
Entrez	10085
Expiration Date	Please enquire.



Representative
Standard
Curve