

Product Datasheet

Human EFEMP1 ELISA Kit (orb385291)

Description Sandwich ELISA for quantitative detection of human Fibulin 3 concentrations.

Reactivity Human

Range 312 pg/ml - 20000 pg/ml

Tested Applications ELISA

Target EFEMP1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q12805**

Sensitivity 10 pg/ml

Entrez **2202**

Expiration Date Please enquire.



Representative
Standard
Curve