


Product Datasheet

Human Hepcidin-25 ELISA Kit (orb385300)

Description	Sandwich ELISA for quantitative detection of human Hepcidin concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	HAMP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P81172
Sensitivity	4 pg/ml
Entrez	57817
Expiration Date	Please enquire.



Representative
Standard
Curve