

Product Datasheet

Human ICAM5 ELISA Kit (orb385303)

Description Sandwich ELISA for quantitative detection of human ICAM-5 concentrations.

Reactivity Human

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target ICAM5

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q9UMF0**

Sensitivity 15 pg/ml

Entrez **7087**

Expiration Date Please enquire.



Representative
Standard
Curve