


Product Datasheet

Human IL17RA ELISA Kit (orb385307)

Description	Sandwich ELISA for quantitative detection of human IL-17RA concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL17RA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q96F46
Sensitivity	10 pg/ml
Entrez	23765
Expiration Date	Please enquire.



Representative
Standard
Curve