

Product Datasheet

Human KLK12 ELISA Kit (orb385310)

Description Sandwich ELISA for quantitative detection of human Kallikrein 12 concentrations.

Reactivity Human

Range 156 pg/ml - 10000 pg/ml

Tested Applications ELISA

Target KLK12

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates

Uniprot ID [Q9UKR0](#)

Sensitivity 10 pg/ml

Entrez [43849](#)

Expiration Date Please enquire.



Representative
Standard
Curve