


Product Datasheet

Human LBP ELISA Kit (orb385312)

Description	Sandwich ELISA for quantitative detection of human LBP concentrations.
Reactivity	Human
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	LBP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P18428
Sensitivity	400 pg/ml
Entrez	3929
Expiration Date	Please enquire.



Representative
Standard
Curve