

Product Datasheet

Human LDLR ELISA Kit (orb385313)

Description Sandwich ELISA for quantitative detection of human LDL Receptor concentrations.

Reactivity Human

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target LDLR

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma, Urine

Uniprot ID **P01130**

Sensitivity 12 pg/ml

Entrez **3949**

Expiration Date Please enquire.



Representative
Standard
Curve