


Product Datasheet

Human NBL1 ELISA Kit (orb385320)

Description	Sandwich ELISA for quantitative detection of human NBL1 concentrations.
Reactivity	Human
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	NBL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P41271
Sensitivity	10 pg/ml
Entrez	4681, 100532736
Expiration Date	Please enquire.



Representative
Standard
Curve