

Product Datasheet

Human NECTIN4 ELISA Kit (orb385321)

Description Sandwich ELISA for quantitative detection of human Nectin 4 concentrations.

Reactivity Human

Range 31.2 pg/ml - 2000 pg/ml

Tested Applications ELISA

Target NECTIN4

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates, Urine

Uniprot ID **Q96NY8**

Sensitivity 10 pg/ml

Entrez **81607**

Expiration Date Please enquire.



Representative
Standard
Curve