

Product Datasheet

Human RLN1 ELISA Kit (orb385331)

Description Sandwich ELISA for quantitative detection of human Relaxin 1 concentrations.

Reactivity Human

Range 78 pg/ml - 5000 pg/ml

Tested Applications ELISA

Target RLN1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P04808**

Sensitivity 10 pg/ml

Entrez **6013**

Expiration Date Please enquire.



Representative
Standard
Curve