

Product Datasheet

Human SAA1 ELISA Kit (orb385338)

Description Sandwich ELISA for quantitative detection of human SAA1 concentrations.

Reactivity Human

Range 1560 pg/ml - 100000 pg/ml

Tested Applications ELISA

Target SAA1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma, Urine

Uniprot ID **P0DJ18**

Sensitivity 150 pg/ml

Entrez **6288**

Expiration Date Please enquire.



Representative
Standard
Curve