

Product Datasheet

Human SHBG ELISA Kit (orb385342)

Description Sandwich ELISA for quantitative detection of human SHBG concentrations.

Reactivity Human

Range 0.031 nmol/L - 2 nmol/L

Tested Applications ELISA

Target SHBG

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **P04278**

Sensitivity 0.001 nmol/L

Entrez **6462**

Expiration Date Please enquire.



Representative
Standard
Curve