

Product Datasheet

Human SPHK1 ELISA Kit (orb385344)

Description Sandwich ELISA for quantitative detection of human SPHK1 concentrations.

Reactivity Human

Range 312 pg/ml - 20000 pg/ml

Tested Applications ELISA

Target SPHK1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID [Q9NYA1](#)

Sensitivity 10 pg/ml

Entrez [8877](#)

Expiration Date Please enquire.



Representative
Standard
Curve