


Product Datasheet

Human SPINT2 ELISA Kit (orb385346)

Description	Sandwich ELISA for quantitative detection of human HAI-2 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	SPINT2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	043291
Sensitivity	10 pg/ml
Entrez	10653
Expiration Date	Please enquire.



Representative
Standard
Curve