


## Product Datasheet

### Human TNFRSF6B ELISA Kit (orb385350)

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sandwich ELISA for quantitative detection of human DCR3 concentrations.                       |
| <b>Reactivity</b>          | Human                                                                                         |
| <b>Range</b>               | 156 pg/ml - 10000 pg/ml                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                         |
| <b>Target</b>              | TNFRSF6B                                                                                      |
| <b>Storage</b>             | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                         |
| <b>Application notes</b>   | Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine            |
| <b>Sample Types</b>        | Cell Culture Supernates, Serum, Plasma                                                        |
| <b>Uniprot ID</b>          | <b>O95407</b>                                                                                 |
| <b>Sensitivity</b>         | 80 pg/ml                                                                                      |
| <b>Entrez</b>              | <b>8771</b>                                                                                   |
| <b>Expiration Date</b>     | Please enquire.                                                                               |



Representative  
Standard  
Curve