

Product Datasheet

Human TNFSF18 ELISA Kit (orb385352)

Description Sandwich ELISA for quantitative detection of human GITRL concentrations.

Reactivity Human

Range 125 pg/ml - 8000 pg/ml

Tested Applications ELISA

Target TNFSF18

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q9UNG2**

Sensitivity 4 pg/ml

Entrez **8995**

Expiration Date Please enquire.



Representative
Standard
Curve