

Product Datasheet

Mouse ACVRL1 ELISA Kit (orb385358)

Description Sandwich ELISA for quantitative detection of mouse ALK1 concentrations.

Reactivity Mouse

Range 46.9 pg/ml - 3000 pg/ml

Tested Applications ELISA

Target ACVRL1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q61288**

Sensitivity 15 pg/ml

Entrez **11482**

Expiration Date Please enquire.



Representative
Standard
Curve