

Product Datasheet

Mouse BCAM ELISA Kit (orb385360)

Description Sandwich ELISA for quantitative detection of mouse CD239 concentrations.

Reactivity Mouse

Range 31.2 pg/ml - 2000 pg/ml

Tested Applications ELISA

Target BCAM

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q9R069**

Sensitivity 15 pg/ml

Entrez **57278**

Expiration Date Please enquire.



Representative
Standard
Curve