

Product Datasheet

Mouse BMP4 ELISA Kit (orb385361)

Description Sandwich ELISA for quantitative detection of mouse BMP4 concentrations.

Reactivity Mouse

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target BMP4

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Bone Tissue, Cell Culture Supernates, Serum

Uniprot ID **P21275**

Sensitivity 10 pg/ml

Entrez **12159**

Expiration Date Please enquire.



Representative
Standard
Curve