


Product Datasheet

Mouse PDL1 ELISA Kit (orb385365)

Description	Mouse PDL1 ELISA Kit
Reactivity	Mouse
Range	125 pg/ml - 8000 pg/ml
Tested Applications	ELISA
Target	CD274
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9EP73
Sensitivity	60 pg/ml
Entrez	60533
Expiration Date	Please enquire.



Representative
Standard
Curve