


Product Datasheet

Mouse CFH ELISA Kit (orb385368)

Description	Sandwich ELISA for quantitative detection of mouse Complement Factor H concentrations.
Reactivity	Mouse
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	CFH
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P06909
Sensitivity	20 pg/ml
Entrez	12628
Expiration Date	Please enquire.



Representative
Standard
Curve