

Product Datasheet

Mouse CXCL3 ELISA Kit (orb385370)

Description Sandwich ELISA for quantitative detection of mouse GRO gamma concentrations.

Reactivity Mouse

Range 15.6 pg/ml - 1000 pg/ml

Tested Applications ELISA

Target CXCL3

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q6W5C0**

Sensitivity 10 pg/ml

Entrez **330122**

Expiration Date Please enquire.



Representative
Standard
Curve