


Product Datasheet

Mouse FABP4 ELISA Kit (orb385375)

Description	Sandwich ELISA for quantitative detection of mouse FABP4 concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	FABP4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	P04117
Sensitivity	15 pg/ml
Entrez	11770
Expiration Date	Please enquire.



Representative
Standard
Curve