


Product Datasheet

Mouse FCGRT ELISA Kit (orb385377)

Description	Sandwich ELISA for quantitative detection of mouse FcRn concentrations.
Reactivity	Mouse
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	FCGRT
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates
Uniprot ID	Q61559
Sensitivity	10 pg/ml
Entrez	14132
Expiration Date	Please enquire.



Representative
Standard
Curve