


Product Datasheet

Mouse LDLR ELISA Kit (orb385386)

Description	Sandwich ELISA for quantitative detection of mouse LDL Receptor concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	LDLR
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	P35951
Sensitivity	12 pg/ml
Entrez	16835
Expiration Date	Please enquire.



Representative
Standard
Curve