


Product Datasheet

Mouse MEP1A ELISA Kit (orb385388)

Description	Sandwich ELISA for quantitative detection of mouse MEP1A concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	MEP1A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P28825
Sensitivity	10 pg/ml
Entrez	17287
Expiration Date	Please enquire.



Representative
Standard
Curve