




Product Datasheet

IRAK4 antibody (orb38539)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	IRAK4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51530
Uniprot ID	Q9NWZ3
NCBI	NP_001138729.1 , NP_001138728.1 , NP_057207.2 , NP_001107654.1 , NP_001138730.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...



Western blot analysis of 293 cell lysate...