

Product Datasheet

Mouse METRNL ELISA Kit (orb385391)

Description Sandwich ELISA for quantitative detection of mouse Meteorin-like concentrations.

Reactivity Mouse

Range 62.5 pg/ml - 4000 pg/ml

Tested Applications ELISA

Target METRNL

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q8VE43**

Sensitivity 12 pg/ml

Entrez **210029**

Expiration Date Please enquire.



Representative
Standard
Curve