


Product Datasheet

Mouse NBL1 ELISA Kit (orb385393)

Description	Sandwich ELISA for quantitative detection of mouse NBL1 concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	NBL1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernates, Serum, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q61477
Sensitivity	10 pg/ml
Entrez	17965
Expiration Date	Please enquire.



Representative
Standard
Curve